

CAEN V1742 Amplitude Calibration

S.Foulkes, S.Los, **P.Murat**, A.Ronzhin, E.Ramberg (Fermilab)


April 4, 2012

Introduction

- DRS4 : two different calibrations: amplitude and time
- here: amplitude calibration only
- DRS4 dynamic range: 1Vpp, symmetric case: +/- 0.5V
- use pulse generator and 5 different voltages: -0.4V, -0.2V, 0V, +0.2V, +0.4V




Calibrating single cell




- width of the line: 3 channels, 1.6 channels w/o generator
- straight line fit through 5 points returns the cell parameters, $\chi^2/N_{DF} \sim 1$


Channels 00-01




Channels 02-03




Channels 04-05




Channels 06-07




Channels 08-09




Channels 10-11




Channels 12-13




Channels 14-15




Channels 16-17




Channels 18-19




Channels 20-21




Channels 22-23









Channels 24-25



Channels 26-27




Channels 28-29









● channels 28-29

Channels 30-31



Summary

Calibration Summary



- no bad channels!
- start using these constants offline
- in parallel: upload the constants into the board - Steve has information from CAEN.