

Studies of the 32-channel TOF-PET detector prototype


Correlated Effects in the Test System Readout

S.Foulkes, S.Los, **P.Murat**, A.Ronzhin, E.Ramberg (Fermilab)

April 5, 2012



Setup

- use new LYSO crystal arrays with 4x4 Hamamatsu 11064 MPPC arrays
- ^{22}Na source (not shown)
- trigger on coincidences of the channels 13 and 21
- $V(\text{MPPC}) = -77\text{V}$, the same on both arrays
- currents very similar: 480mA and 474 mA



Introduction

- V1742 calibrated, “per-cell” pedestal subtraction works
- resulting noise RMS < 1mV
- pulse charge $Q = \Sigma V_i / 10$
- expect small charge in empty channels, this, however, depends on the pedestal subtraction
- pedestal shifts \rightarrow large $|Q|$
- below: a trivial event display, showing charges in the two 4x4 crystal arrays
- left 4x4 channels: the first array, right 4x4 channels: the second array




- pulses in the two highest channels


An Example of a “Normal” Event

- so far: calculate and subtract average pedestals
- shown event: result of subtraction is consistent with zero


Event Q



v_014




v_020





Correlated shifts of the pedestals

- observe events with the zero voltage level shifted wrt zero
- the shift is correlated across the channels



- to account for this effect, need to fit pedestals on event-by-event basis in each channel using cells prior to the pulse



- pulses in several channels

Summary

- amplitude calibrations look good, remaining noise RMS < 1 mV
- pulses are good, S/N > 100
- observe pedestal shifts correlated across the channels
- need to fit pedestals on event-by-event basis
- currently, T0 ~ 500 channels, need to reduce the trigger delay
- want to leave ~ 100 channels (20ns) for fitting the pedestal
- with the current signals, it is sufficient to digitize and read out 400 channels = 80ns